

Material	Dk	Wetting Angle	Hardness	Handling Tint	Surface	Available MNL Designs
Optimum Comfort	65	6	78	Green/Blue	Not surface treated as naturally wettable	All
Paragon HDS (Futura)	58	15	84	Green/Blue	Not surface treated as naturally wettable	All
Menicon EX	85	14	75	Blue	Plasma Treated to help with comfort and wettability but only last 10 days approx	All
Menicon Z	189	22	83	Blue	Plasma Treated to help with comfort and wettability but only last 10 days approx	All
Boston ES	18	52	85	Blue	Not surface treated as naturally wettable	All Omni/Exact/FlexiCurve/Comfort Symbio and Progressive Near only
Boston XO	100	49	79	Ice Blue	Not surface treated as naturally wettable	All Comfort/Omni/Exact/Flexicurve/Rose K2 / Rose K2 IC and Rose K2 PG
Boston XO	100	49	79	Violet	Not surface treated as naturally wettable	All Comfort/Omni/Exact and Flexicurve Rose K2 / Rose K2 IC/ Rose K2 NC and Rose K2 PG